

***BILI-THERAPY** Pad Type*

LED Phototherapy Unit




CE 0123




BILI-THERAPY Pad Type delivers effective phototherapy for the treatment of indirect hyperbilirubinemia from the back of the newborn.

•Relative irradiance



Effective wavelength from the light source of blue LED will encourage newborn's early discharge from the hospital.



BILI-THERAPY Pad Type can be used in a bassinet next to the mother's bed.



BILI-THERAPY Pad Type along with BILI-THERAPY Spot Type can effectively provide double phototherapy treatment without interfering with visualisation of the newborn.

Specifications

BILI-THERAPY Pad Type	
Atom's code	60370
Power requirement	Customer specified
Power consumption	68VA
Frequency	50/60Hz
Type of protection	Type BF applied part
Degree of protection	Class I
Classification	IIa
Dimensions	150(W)x220(D)x160(H) mm
Weight	Approx. 2.1Kg
Accessories	Power cord (AC120V) x 1, Power cord (AC230V) x 1, Operation manual x 1
Light source	Blue LEDs
Irradiance level	High : 53μW/cm ² /nm (3.2mW/cm ²)±25% Medium : 40μW/cm ² /nm (2.4mW/cm ²)±25% Low : 26.5μW/cm ² /nm (1.6mW/cm ²)±25% obtained by the BiliBlanket® Meter at the irradiation distance of 30cm

Options

Atom's code	Description
60371	Pad (L) For BILI-THERAPY Pad Type
60372	Pad (S) For BILI-THERAPY Pad Type
60373	Pad Cover (L) For BILI-THERAPY Pad Type (50 pcs/box)
60374	Pad Cover (S) For BILI-THERAPY Pad Type (50 pcs/box)
60375	BILI-Holder For BILI-THERAPY Pad Type
60217	F-Rail Attachment
60218	Pole Attachment
60376	BILI-Stand For BILI-THERAPY Pad Type
68448	Phototherapy Analyzer II

Manufacturer: ATOM MEDICAL CORPORATION

18-15, Hongo 3-chome, Bunkyo-ku, Tokyo, 113-0033 JAPAN
<http://www.atomed.co.jp>

Exporter: ATOM MEDICAL INTERNATIONAL, INC.

SWT Bldg., 15-9 Hongo 3-chome, Bunkyo-ku, Tokyo, 113-0033 JAPAN
 Phone: + 81 3 3815 2941 Fax: + 81 3 3812 9670/4080

BPL MEDICAL TECHNOLOGIES PRIVATE LIMITED

11th KM, Bannerghatta Road, Arakere, Bangalore - 560076
 CIN U33110KA2012PTC067282
 Phone +91 80 26484388/ 2648 4350
 Toll Free 1800-425-2355
 Website www.bplmedicaltechnologies.com
 For Enquiries: sales.medical@bpl.in
 Copyright©2017 BPL Medical Technologies Private Limited.
 All rights reserved.