

Experience the Flexibility

BPL120/80_{B1}**Automatic Blood Pressure Monitor***Happier Living Everyday*

Pressure Check Made Easy

**Automatic Pressurization**

High-tech measurement technology to monitor blood pressure from the arm – all with the push of a button

**Large LCD Display**

Extra large LCD display of 84.1 X 55.1 mm provides excellent display clarity and enables easy reading

**Irregular Heartbeat Detection**

Automatically alerts the user when heart rhythm is either 25% faster or slower than the average rhythm

**WHO BP Classification Indicator**

Automatically indicates blood pressure classification after the completion of each measurement

**120 Memories in 4 Groups**

Automatic memory function records the last 120 measurements in 4 groups for efficient and simplified tracking

**Low Battery Detection**

Visually indicates when the battery has been discharged below a specific level

CERTIFIED ISO 13485:2003, ISO 9001:2008 COMPANY

BPL Medical Technologies Private Limited
11th KM, Bannerghatta Road,
Arakere, Bangalore - 560076
CIN: U33110KA2012PTC067282
Ph: +91 80 26484388/ 2648 4350
Toll Free: 1800-4252355
Website: www.bplmedicaltechnologies.com
For Enquiries: sales.medical@bpl.in



A BPL Medical Technologies Magazine

