



Experience the Trace Quality



Happier Living Everyday

THREE CHANNEL ECG RECORDER **CARDIART 7108**



Three Channel ECG Recording

Simultaneous 3 (user selectable – 3/2/1 channel) channel ECG recording



Last Memorised ECG Recording

Ensures the storage of last recording in the memory



Automatic and Manual Recording Modes

Manual mode enables the acquisition and printing of ECG recording in real time. Simultaneous acquisition is possible with automatic mode



Power Supply

Comes with dual power supply – mains and rechargeable internal battery



Built-in Rechargeable Battery

Enables mains independent operations and provides battery use after charging



Thermal Printing

Allows the printing of ECG on 60 mm paper using a high-resolution thermal print



Ergonomic Design

The compact, light weight and portable design ensures ease of use



Printing of Rhythm Lead Trace

Printing of user-selectable rhythm lead trace up to 100 seconds



Optional

Stored and real-time ECG transfer to PC through USB

CERTIFIED ISO 13485:2003, ISO 9001:2008 COMPANY

BPL Medical Technologies Private Limited
11th KM, Bannerghatta Road,
Arakere, Bangalore - 560076
CIN: U33110KA2012PTC067282
Ph: +91 80 26484388/ 2648 4350
Toll Free: 1800-4252355
Website: www.bplmedicaltechnologies.com
For Enquiries: sales.medical@bpl.in



A BPL Medical Technologies
Magazine



Follow Us on

