

BPL H-RAD₃₂

BPL Medical
Technologies

The Comprehensive High Frequency Fixed X-Ray Imaging solution



www.bplmedicaltechnologies.com



Happier Living Everyday



Follow us on



A BPL Medical Technologies
Magazine

Powered by



Scan the product image
above with BPL AR App to
view the product video of
BPL H-Rad 32



Download BPL-AR App to see this catalogue come to life

Download BPL Promise App to experience the latest in
medical technology

BPL H-RAD₃₂

BPL H-Rad 32 is a 32 kW High Frequency Fixed X-Ray solution that combines excellent image quality and comprehensive User Interface. It is suitable for wide Clinical applications and is ideal for Govt Hospitals, Medical institutions, Super / Multi Speciality Hospitals, Corporate Hospitals, Diagnostic Centres etc.,


The Comprehensive High Frequency Fixed X-Ray Imaging solution

- 1 User Friendly Design
- 2 Superior Image Quality
- 3 Maximum Value
- 4 Minimum Downtime



Powered By: High Frequency 32 kW X-Ray Generator

BPL H Rad 32 is powered by High Frequency 32 kW X-Ray Generator and is equipped with an Advanced User friendly Operator Console with soft touch control of all radiographic parameters.



Approved by



- User Configurable Anatomical Programming with 216 Program options
- Self-Diagnostic Program
- 40 to 125 kVp in steps of 1 kVp only
- Density Control
- Automatic Exposure Control (Optional)



Easy to use Counter Balanced Tube Head




Automatic Calibration helps in high consistency in exposures and image quality by reproducing Exposure parameters using the feedback loop (Optional).




ERROR

Self-Diagnostic Program to minimize downtime



Anatomical Programming for precise selection of kV & mAs based on patient's anatomy (216 Program options)




Remote Diagnostic & Service Tool (Optional)



Automatic Exposure Control to optimize X-Ray Dose to Patient & Operator (Optional)

A Global Medical Technology Company



Service is Powered by



Service Helpline
1800-425-2355

CERTIFIED ISO 13485:2003, ISO 9001:2008 COMPANY

BPL Medical Technologies Private Limited
11th KM, Bannerghatta Road,
Arakere, Bangalore - 560076
CIN: U33110KA2012PTC067282
Ph: +91 80 26484388/ 2648 4350
Toll Free: 1800-4252355
Website: www.bplmedicaltechnologies.com
For Enquiries: sales.medical@bpl.in



ALWAYS USE GENUINE BPL
CONSUMABLES AND
ACCESSORIES FOR
BETTER PERFORMANCE



Happier Living Everyday