



Experience the Compactness



Happier Living Everyday

Experience the freedom of easy and continuous acquisition of high-quality ECG readings with utmost convenience



Compact and Lightweight

68mmX53mmX16mm. 42 grams (excluding battery)



Fully Programmable, High Fidelity, High Quality Performance

Sample rate up to 1024 samples per second/ no data compression



12-Lead/ 3 Channel Two Work Modes

Two work modes are available:
10 Wire Patient Cable (12 Leads)
5 Wire Patient Cable/ 7 Wire Patient Cable (3 Lead)



Lower Power Consumption

Powered by an AAA size alkaline/ NiMH battery > 48 hours



Removable SD Card

2GB SD card records up to 48 hours of the 12-lead ECG at 1024 SPS



LCD Display

Provides on demand ECG view/ controls programming instructions/ visual confirmation of patient hook up



Pacemaker Detection

Can detect the small, narrow pulses from the pacemaker



Two Optional Modes of ECG Data transfer

Connected directly with PC through USB connector embedded in the recorder/ removable SD card

CERTIFIED ISO 13485:2003, ISO 9001:2008 COMPANY

BPL Medical Technologies Private Limited
11th KM, Bannerghatta Road,
Arakere, Bangalore - 560076
CIN: U33110KA2012PTC067282
Ph: +91 80 26484388/ 2648 4350
Toll Free: 1800-4252355
Website: www.bplmedicaltechnologies.com
For Enquiries: sales.medical@bpl.in



A BPL Medical Technologies Magazine



Follow Us on

