

The Reliable Mobile X-Ray Imaging Solution



www.bplmedicaltechnologies.com



Happier Living Everyday



Follow us on



A BPL Medical Technologies
Magazine

Powered by



Scan the product image
above with BPL AR App to
view the product video of
BPL M Rad 100



Download BPL-AR App to see this catalogue come to life

Download BPL Promise App to experience the latest in
medical technology

BPL M RAD™ 100

BPL M Rad 100 is a 100 mA, 100 PPS Line Frequency Mobile X-ray solution that combines Excellent Manoeuvrability & User Interface and is ideal for patient wards, intensive care, operating room etc., where transporting the Patient to the Radiology Department is difficult



Approved by




The Reliable Mobile X-Ray Imaging solution

- 1 User Friendly Design
- 2 Superior Image Quality
- 3 Maximum Value
- 4 Minimum Downtime

Powered By:



BPL M Rad 100 is powered by Advanced User friendly Operator Console with soft touch control of all radiographic parameters. The equipment uses State of the Art Microcontroller to control various functions and features.



- Alphanumeric LCD Display (20X4 Characters)
- User Configurable Anatomical Programming
- 155 Anatomical Programs
- Pre-Set Key for instant selection
- Self-Diagnostic Program
- X-Ray Counter with password protection
- Memory to retain last selected factors
- 40 to 100 kVp in steps of 2 kVp only



Counter Balanced Tube Head with Vertical, Horizontal, Rotational, Axial & Swivel movements




Durable wheels & Efficient Foot Brake mechanism




Storage space for Accessories, Notepads etc.



Independent Voltmeter for monitoring input line Voltage & independent Ammeter for indication of Tube Current




Anatomical Programming for precise selection of kV & mAs based on patient's anatomy



Self-Diagnostic Program to minimize downtime

A Global Medical Technology Company



Service is Powered by



Service Helpline
1800-425-2355

CERTIFIED ISO 13485:2003, ISO 9001:2008 COMPANY

BPL Medical Technologies Private Limited
11th KM, Bannerghatta Road,
Arakere, Bangalore - 560076
CIN: U33110KA2012PTC067282
Ph: +91 80 26484388/ 2648 4350
Toll Free: 1800-4252355
Website: www.bplmedicaltechnologies.com
For Enquiries: sales.medical@bpl.in



ALWAYS USE GENUINE BPL
CONSUMABLES AND
ACCESSORIES FOR
BETTER PERFORMANCE



Happier Living Everyday