



Experience the Quality Care

*Happier Living Everyday*

MATERNAL MONITOR

PRIMA-MM 9855



Monitor Screen

12.1" high resolution colour TFT display with tiltable screen up to 90°



Standard Configuration

FHR, TOCO, Foetal movement, ECG, SpO₂, NiBP, RESP and TEMP



Portability

The device is compact, portable and lightweight



Ultrasound Probe

The probe with high sensitivity enables user to monitors foetus > 12 weeks and provide versatility in clinical settings



Data Display & Storage

The device can record and display data for 25 hours with dynamic data save capability



Thermal Printer

Thermal printer facilitates easy availability of printed results



Safety

The low frequency ultrasound waves are safe on foetus



AFM

Automatic Foetal Movement detection with event marker



Li-ion Battery

The built-in rechargeable Li-ion battery ensures uninterrupted use of the device



Optional

Twin FHR Monitoring & Central Monitoring System Software

CERTIFIED ISO 13485:2003, ISO 9001:2008 COMPANY

BPL Medical Technologies Private Limited
11th KM, Bannerghatta Road,
Arakere, Bangalore - 560076
CIN: U33110KA2012PTC067282
Ph: +91 80 26484388/ 2648 4350
Toll Free: 1800-4252355
Website: www.bplmedicaltechnologies.com
For Enquiries: sales.medical@bpl.in



A BPL Medical Technologies
Magazine



Follow Us on

