



Experience the Effectiveness

BIPHASIC DEFIBRILLATOR
RELIFE 700



Happier Living Everyday



Saving lives through enhanced technology with improved safety and effectiveness



Biphasic Technology

The current-controlled biphasic waveform lowers the electrical threshold for successful defibrillation with minimal tissue damage and the truncated exponential waveform is adapted to patient's impedance



Simple and Intuitive

The machine is user-friendly with simple and intuitive interface



Portable and Lightweight

The machine is handy and easy to carry



Remote Control through Paddles

Remote charging and discharging through paddles facilitates easy charging for convenient use



Integrated Paediatric Paddles

Designed for safe use on paediatric patients enabling suitable doses of energy



Flexible Power Supply

Operates from the mains(AC), from a vehicle battery(DC) and with its internal battery



Autotest

Enables autotest when switched on and during operations



Optional

AED, Pacemaker, MASIMO SpO₂, Recorder and Data storage

CERTIFIED ISO 13485:2003, ISO 9001:2008 COMPANY

BPL Medical Technologies Private Limited
11th KM, Bannerghatta Road,
Arakere, Bangalore - 560076
CIN: U33110KA2012PTC067282
Ph: +91 80 26484388/ 2648 4350
Toll Free: 1800-4252355
Website: www.bplmedicaltechnologies.com
For Enquiries: sales.medical@bpl.in



A BPL Medical Technologies Magazine



Follow Us on

